



F Antibody

Product Code	CSB-PA318261LA01CCQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12569
Immunogen	Recombinant Canine distemper virus Fusion glycoprotein F0 protein (136-608AA)
Raised In	Rabbit
Species Reactivity	Canine distemper virus
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fusion glycoprotein F0 [Cleaved into: Fusion glycoprotein F2; Fusion glycoprotein F1], F
Immunogen Species	Canine distemper virus (strain Onderstepoort) (CDV)
Research Area	Others
Target Names	F