



# mgsA Antibody, Biotin conjugated

Product Code	CSB-PA507015HD01BQK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C1EN30
Immunogen	Recombinant Bacillus cereus Methylglyoxal synthase protein (1-131AA)
Raised In	Rabbit
Species Reactivity	Bacillus cereus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methylglyoxal synthase (MGS) (EC 4.2.3.3), mgsA
Immunogen Species	Bacillus cereus (strain 03BB102)
Research Area	Others
Target Names	mgsA