



hly Antibody, Biotin conjugated

Product Code	CSB-PA639324HD01FLF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2G1X0
Immunogen	Recombinant Staphylococcus aureus Alpha-hemolysin protein (27-319AA)
Raised In	Rabbit
Species Reactivity	Staphylococcus aureus
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-hemolysin (Alpha-HL) (Alpha-toxin), hly, hla
Immunogen Species	Staphylococcus aureus (strain NCTC 8325)
Research Area	Others
Target Names	hly
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.