



VTI1B Antibody

Product Code	CSB-PA887014LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UEU0
Immunogen	Recombinant Human Vesicle transport through interaction with t-SNAREs homolog 1B protein (2-208AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Vesicle transport through interaction with t-SNAREs homolog 1B (Vesicle transport v-SNARE protein Vti1-like 1) (Vti1-rp1), VTI1B, VT11 VTI1L VTI1L1 VTI2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	VTI1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.