



RORC Antibody, FITC conjugated

Product Code	CSB-PA020071LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51449
Immunogen	Recombinant Human Nuclear receptor ROR-gamma protein (1-518AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Nuclear receptor ROR-gamma (Nuclear receptor RZR-gamma) (Nuclear receptor subfamily 1 group F member 3) (RAR-related orphan receptor C) (Retinoid-related orphan receptor-gamma), RORC, NR1F3 RORG RZRG
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	RORC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.