



Arachin Ahy-3 Antibody

Product Code	CSB-PA737391LA01ANE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q647H2
Immunogen	Recombinant Peanut Arachin Ahy-3 protein (299-478AA)
Raised In	Rabbit
Species Reactivity	Peanut
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Arachin Ahy-3 [Cleaved into: Arachin Ahy-3 chain alpha; Arachin Ahy-3 chain beta]
Immunogen Species	Peanut
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.