



# Sort1 Antibody

Product Code	CSB-PA022412LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O54861
Immunogen	Recombinant Rat Sortilin protein (610-754AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sortilin (Glycoprotein 110) (Gp110) (Neurotensin receptor 3) (NTR3), Sort1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Sort1