



Sort1 Antibody, HRP conjugated

Product Code	CSB-PA022412LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O54861
Immunogen	Recombinant Rat Sortilin protein (610-754AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sortilin (Glycoprotein 110) (Gp110) (Neurotensin receptor 3) (NTR3), Sort1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Neuroscience
Target Names	Sort1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.