



IL10 Antibody, Biotin conjugated

Product Code	CSB-PA011580LD01BO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43480
Immunogen	Recombinant Bovine Interleukin-10 protein (19-178AA)
Raised In	Rabbit
Species Reactivity	Bovine
Tested Applications	ELISA
Relevance	Inhibits the synthesis of a number of cytokines, including IFN-gamma, IL-2, IL-3, TNF and GM-CSF produced by activated macrophages and by helper T-cells.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-10 (IL-10) (Cytokine synthesis inhibitory factor) (CSIF), IL10
Immunogen Species	Bos taurus (Bovine)
Research Area	Others
Target Names	IL10