



OGT Antibody, FITC conjugated

Product Code	CSB-PA016315LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15294
Immunogen	Recombinant Human UDP-N-acetylglucosamine--peptide N-acetylglucosaminyltransferase 110 kDa subunit protein (606-1022AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UDP-N-acetylglucosamine--peptide N-acetylglucosaminyltransferase 110 kDa subunit (EC 2.4.1.255) (O-GlcNAc transferase subunit p110) (O-linked N-acetylglucosamine transferase 110 kDa subunit) (OGT), OGT
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	OGT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.