



TYR Antibody, HRP conjugated

Product Code	CSB-PA025394LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14679
Immunogen	Recombinant Human Tyrosinase protein (19-377AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosinase (EC 1.14.18.1) (LB24-AB) (Monophenol monooxygenase) (SK29-AB) (Tumor rejection antigen AB), TYR
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	TYR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.