



Selenot Antibody, HRP conjugated

Product Code	CSB-PA350436LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62342
Immunogen	Recombinant Mouse Thioredoxin reductase-like selenoprotein T protein (20-195AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thioredoxin reductase-like selenoprotein T (SeIT) (EC 1.8.1.9), Selenot, Selt
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Selenot
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.