



ompA Antibody, Biotin conjugated

Product Code	CSB-PA328425LD01DSA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23732
Immunogen	Recombinant Chlamydia trachomatis Major outer membrane porin, serovar A protein (23-396AA)
Raised In	Rabbit
Species Reactivity	Chlamydia trachomatis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane porin, serovar A (MOMP), ompA, omp1A
Immunogen Species	Chlamydia trachomatis
Research Area	Others
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.