



FCER2 Antibody, HRP conjugated

Product Code	CSB-PA02219B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06734
Immunogen	Recombinant Human Low affinity immunoglobulin epsilon Fc receptor protein (48-321AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin epsilon Fc receptor (BLAST-2) (C-type lectin domain family 4 member J) (Fc-epsilon-R1I) (Immunoglobulin E-binding factor) (Lymphocyte IgE receptor) (CD antigen CD23) [Cleaved into: Low affinity immunoglobulin epsilon Fc receptor membrane-bound form; Low affinity immunoglobulin epsilon Fc receptor soluble form], FCER2, CD23A CLEC4J FCE2 IGEBF
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCER2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.