



ompA Antibody, Biotin conjugated

Product Code	CSB-PA333258LD01DRZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P27455
Immunogen	Recombinant Chlamydomonas reinhardtii Major outer membrane porin protein (24-389AA)
Raised In	Rabbit
Species Reactivity	Chlamydomonas reinhardtii
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Major outer membrane porin (MOMP), ompA, omp1
Immunogen Species	Chlamydia pneumoniae (Chlamydomonas reinhardtii)
Research Area	Others
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.