



NPR2 Antibody, Biotin conjugated

Product Code	CSB-PA016024LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P20594
Immunogen	Recombinant Human Atrial natriuretic peptide receptor 2 protein (23-458AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Atrial natriuretic peptide receptor 2 (EC 4.6.1.2) (Atrial natriuretic peptide receptor type B) (ANP-B) (ANPR-B) (NPR-B) (Guanylate cyclase B) (GC-B), NPR2, ANPRB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NPR2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.