



bla Antibody, HRP conjugated

Product Code	CSB-PA352353YB01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62593
Immunogen	Recombinant Escherichia coli Beta-lactamase TEM protein (24-286AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Beta-lactamase TEM (EC 3.5.2.6) (IRT-4) (Penicillinase) (TEM-1) (TEM-16/CAZ-7) (TEM-2) (TEM-24/CAZ-6) (TEM-3) (TEM-4) (TEM-5) (TEM-6) (TEM-8/CAZ-2), bla; blaT-3; blaT-4; blaT-5; blaT-6,
Immunogen Species	Escherichia coli
Research Area	Others
Target Names	bla
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.