



# bla Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA352353YD01ENL
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P62593
<b>Immunogen</b>	Recombinant Escherichia coli Beta-lactamase TEM protein (24-286AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Beta-lactamase TEM (EC 3.5.2.6) (IRT-4) (Penicillinase) (TEM-1) (TEM-16/CAZ-7) (TEM-2) (TEM-24/CAZ-6) (TEM-3) (TEM-4) (TEM-5) (TEM-6) (TEM-8/CAZ-2), bla; blaT-3; blaT-4; blaT-5; blaT-6,
<b>Immunogen Species</b>	Escherichia coli
<b>Research Area</b>	Others
<b>Target Names</b>	bla
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.