



dsdA Antibody, Biotin conjugated

Product Code	CSB-PA472744HD01SWQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B4TA53
Immunogen	Recombinant Salmonella heidelberg D-serine dehydratase protein (1-440AA)
Raised In	Rabbit
Species Reactivity	Salmonella heidelberg
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	D-serine dehydratase (EC 4.3.1.18) (D-serine deaminase) (DSD), dsdA
Immunogen Species	Salmonella heidelberg (strain SL476)
Research Area	Others
Target Names	dsdA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.