



ORF3 Antibody

Product Code	CSB-PA527126LA01HVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O90299
Immunogen	Recombinant Hepatitis E virus genotype 1 Protein ORF3 protein (1-114AA)
Raised In	Rabbit
Species Reactivity	Hepatitis E virus genotype 1
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ORF3 (pORF3), ORF3
Immunogen Species	Hepatitis E virus genotype 1
Research Area	Others
Target Names	ORF3