




PADI2 Antibody

Product Code	CSB-PA896493LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2J8
Immunogen	Recombinant Human Protein-arginine deiminase type-2 protein (1-665AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-arginine deiminase type-2 (EC 3.5.3.15) (PAD-H19) (Peptidylarginine deiminase II) (Protein-arginine deiminase type II), PADI2, KIAA0994 PAD2 PD12
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PADI2

Image



Western blot

All lanes: Protein-arginine deiminase type-2 antibody at 16µg/ml

Lane 1: PC-3 whole cell lysate

Lane 2: Hela whole cell lysate

Lane 3: Mouse brain tissue

Lane 4: MDA-MB-231 whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 76, 50 kDa

Observed band size: 76 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.