



PADI2 Antibody, Biotin conjugated

Product Code	CSB-PA896493LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2J8
Immunogen	Recombinant Human Protein-arginine deiminase type-2 protein (1-665AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein-arginine deiminase type-2 (EC 3.5.3.15) (PAD-H19) (Peptidylarginine deiminase II) (Protein-arginine deiminase type II), PADI2, KIAA0994 PAD2 PD12
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PADI2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.