



GPT Antibody, FITC conjugated

Product Code	CSB-PA009864LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24298
Immunogen	Recombinant Human Alanine aminotransferase 1 protein (2-496AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alanine aminotransferase 1 (ALT1) (EC 2.6.1.2) (Glutamate pyruvate transaminase 1) (GPT 1) (Glutamic--alanine transaminase 1) (Glutamic--pyruvic transaminase 1), GPT, AAT1 GPT1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GPT
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.