




TGFB2 Antibody

Product Code	CSB-PA07319A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P61812
Immunogen	Recombinant Human Transforming growth factor beta-2 proprotein protein (304-413AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transforming growth factor beta-2 (TGF-beta-2) (BSC-1 cell growth inhibitor) (Cetermin) (Glioblastoma-derived T-cell suppressor factor) (G-TSF) (Polyergin) [Cleaved into: Latency-associated peptide (LAP)], TGFB2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TGFB2

Image



Western blot

All lanes: TGFB2 antibody at 3µg/ml

Lane 1: Mouse heart tissue


Lane 2: Mouse skeletal muscle tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 48, 51 kDa

Observed band size: 48 kDa



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA07319A0Rb at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.