



p46 Antibody, Biotin conjugated

Product Code	CSB-PA314621LD01MLQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C0J8
Immunogen	Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen protein (28-416AA)
Raised In	Rabbit
Species Reactivity	Mycoplasma hyopneumoniae
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	46 kDa surface antigen (p46), p46
Immunogen Species	Mycoplasma hyopneumoniae (strain J / ATCC 25934 / NCTC 10110)
Research Area	Others
Target Names	p46
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.