



# p46 Antibody, Biotin conjugated

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | CSB-PA314621YD01MLQ  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| <b>Uniprot No.</b>         | P0C0J8   |
| <b>Immunogen</b>           | Recombinant Mycoplasma hyopneumoniae 46 kDa surface antigen protein (28-416AA)   |
| <b>Raised In</b>           | Rabbit   |
| <b>Species Reactivity</b>  | Mycoplasma hyopneumoniae   |
| <b>Tested Applications</b> | ELISA  |
| <b>Form</b>                | Liquid   |
| <b>Conjugate</b>           | Biotin   |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| <b>Purification Method</b> | >95%, Protein G purified   |
| <b>Isotype</b>             | IgG  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Alias</b>               | 46 kDa surface antigen (p46), p46  |
| <b>Immunogen Species</b>   | Mycoplasma hyopneumoniae (strain J / ATCC 25934 / NCTC 10110)                    |
| <b>Research Area</b>       | Others   |
| <b>Target Names</b>        | p46  |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.      |