



BCL2L1 Antibody, Biotin conjugated

Product Code	CSB-PA002613HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q07817
Immunogen	Recombinant Human Bcl-2-like protein 1 protein (1-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bcl-2-like protein 1 (Bcl2-L-1) (Apoptosis regulator Bcl-X), BCL2L1, BCL2L BCLX
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BCL2L1