





# SNCG Antibody

<b>Product Code</b>	CSB-PA021915HA01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	O76070
<b>Immunogen</b>	Recombinant Human Gamma-synuclein protein (1-127AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human, Mouse
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Gamma-synuclein (Breast cancer-specific gene 1 protein) (Persyn) (Synoretin) (SR), SNCG, BCSG1 PERSYN PRSN
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Neuroscience
<b>Target Names</b>	SNCG


## Image



Western blot  
 All lanes: SNCG antibody at 16µg/ml  
 Lane 1: HT29 whole cell lysate  
 Lane 2: Mouse brain tissue  
 Lane 3: 293T whole cell lysate  
 Secondary  
 Goat polyclonal to rabbit IgG at 1/10000 dilution  
 Predicted band size: 14 kDa  
 Observed band size: 14 kDa



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA021915HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA021915HA01HU at dilution of 1:100