





SNCG Antibody

Product Code	CSB-PA021915HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O76070
Immunogen	Recombinant Human Gamma-synuclein protein (1-127AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Gamma-synuclein (Breast cancer-specific gene 1 protein) (Persyn) (Synoretin) (SR), SNCG, BCSG1 PERSYN PRSN
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNCG


Image



Western blot
 All lanes: SNCG antibody at 16µg/ml
 Lane 1: HT29 whole cell lysate
 Lane 2: Mouse brain tissue
 Lane 3: 293T whole cell lysate
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 14 kDa
 Observed band size: 14 kDa



Immunohistochemistry of paraffin-embedded human ovarian cancer using CSB-PA021915HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA021915HA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.