



Streptavidin Antibody, Biotin conjugated

Product Code	CSB-PA332849YD01SNO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22629
Immunogen	Recombinant Streptomyces avidinii Streptavidin protein (25-183AA)
Raised In	Rabbit
Species Reactivity	Streptomyces avidinii
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Streptavidin
Immunogen Species	Streptomyces avidinii
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.