



IL10 Antibody

Product Code	CSB-PA011580LA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51496
Immunogen	Recombinant Macaca mulatta Interleukin-10 protein (19-178AA)
Raised In	Rabbit
Species Reactivity	Macaca mulatta
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-10 (IL-10) (Cytokine synthesis inhibitory factor) (CSIF), IL10
Immunogen Species	Macaca mulatta (Rhesus macaque)
Research Area	Others
Target Names	IL10