



IL10 Antibody

Product Code	CSB-PA011580LA01MOW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P51496
Immunogen	Recombinant Macaca mulatta Interleukin-10 protein (19-178AA)
Raised In	Rabbit
Species Reactivity	Macaca mulatta
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-10 (IL-10) (Cytokine synthesis inhibitory factor) (CSIF), IL10
Immunogen Species	Macaca mulatta (Rhesus macaque)
Research Area	Others
Target Names	IL10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.