



IL10 Antibody, Biotin conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA011580LD01MOW |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P51496 |
| Immunogen | Recombinant Macaca mulatta Interleukin-10 protein (19-178AA) |
| Raised In | Rabbit |
| Species Reactivity | Macaca mulatta |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Interleukin-10 (IL-10) (Cytokine synthesis inhibitory factor) (CSIF), IL10 |
| Immunogen Species | Macaca mulatta (Rhesus macaque) |
| Research Area | Others |
| Target Names | IL10 |