



# prn Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA321224LC01BUA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P14283
<b>Immunogen</b>	Recombinant Bordetella pertussis Pertactin autotransporter protein (632-910AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Bordetella pertussis
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Pertactin autotransporter (P.93) [Cleaved into: Outer membrane protein P.69; Pertactin translocator], prn, omp69A
<b>Immunogen Species</b>	Bordetella pertussis (strain Tohama I / ATCC BAA-589 / NCTC 13251)
<b>Research Area</b>	Others
<b>Target Names</b>	prn
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.