



Map2k4 Antibody, HRP conjugated

Product Code	CSB-PA013413LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47809
Immunogen	Recombinant Mouse Dual specificity mitogen-activated protein kinase kinase 4 protein (2-397AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dual specificity mitogen-activated protein kinase kinase 4 (MAP kinase kinase 4) (MAPKK 4) (EC 2.7.12.2) (C-JUN N-terminal kinase kinase 1) (JNK kinase 1) (JNKK 1) (JNK-activating kinase 1) (MAPK/ERK kinase 4) (MEK 4) (SAPK/ERK kinase 1) (SEK1), Map2k4, Jnkk1 Mek4 Mkk4 Prkmk4 Sek1 Serk1 Skk1
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Map2k4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.