



# Map2k4 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA013413LB01MO
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P47809
<b>Immunogen</b>	Recombinant Mouse Dual specificity mitogen-activated protein kinase kinase 4 protein (2-397AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Dual specificity mitogen-activated protein kinase kinase 4 (MAP kinase kinase 4) (MAPKK 4) (EC 2.7.12.2) (C-JUN N-terminal kinase kinase 1) (JNK kinase 1) (JNKK 1) (JNK-activating kinase 1) (MAPK/ERK kinase 4) (MEK 4) (SAPK/ERK kinase 1) (SEK1), Map2k4, Jnkk1 Mek4 Mkk4 Prkmk4 Sek1 Serk1 Skk1
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	Map2k4