



Ctss Antibody, HRP conjugated

Product Code	CSB-PA006204LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O70370
Immunogen	Recombinant Mouse Cathepsin S protein (123-340AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cathepsin S (EC 3.4.22.27), Ctss, Cats
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Ctss