



SNCA Antibody, FITC conjugated

Product Code	CSB-PA06997C0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37840
Immunogen	Recombinant Human Alpha-synuclein protein (1-140AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-synuclein (Non-A beta component of AD amyloid) (Non-A4 component of amyloid precursor) (NACP), SNCA, NACP PARK1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SNCA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.