



RMDN3 Antibody, HRP conjugated

Product Code	CSB-PA836307LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96TC7
Immunogen	Recombinant Human Regulator of microtubule dynamics protein 3 protein (1-470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of microtubule dynamics protein 3 (RMD-3) (hRMD-3) (Cerebral protein 10) (Protein FAM82A2) (Protein FAM82C) (Protein tyrosine phosphatase-interacting protein 51) (TCPTP-interacting protein 51), RMDN3, FAM82A2 FAM82C PTPIP51
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RMDN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.