



RMDN3 Antibody, Biotin conjugated

Product Code	CSB-PA836307LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96TC7
Immunogen	Recombinant Human Regulator of microtubule dynamics protein 3 protein (1-470AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Regulator of microtubule dynamics protein 3 (RMD-3) (hRMD-3) (Cerebral protein 10) (Protein FAM82A2) (Protein FAM82C) (Protein tyrosine phosphatase-interacting protein 51) (TCPTP-interacting protein 51), RMDN3, FAM82A2 FAM82C PTPIP51
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RMDN3