



TP53 Antibody, Biotin conjugated

Product Code	CSB-PA07889D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04637
Immunogen	Recombinant Human Cellular tumor antigen p53 protein (2-393AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53), TP53, P53
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	TP53
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.