



JPT1 Antibody, FITC conjugated

Product Code	CSB-PA890732HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UK76
Immunogen	Recombinant Human Jupiter microtubule associated homolog 1 protein (2-154AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Jupiter microtubule associated homolog 1 (Androgen-regulated protein 2) (Hematological and neurological expressed 1 protein) [Cleaved into: Jupiter microtubule associated homolog 1, N-terminally processed], JPT1, ARM2 HN1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	JPT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.