



PNP Antibody, HRP conjugated

Product Code	CSB-PA355527HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P00491
Immunogen	Recombinant Human Purine nucleoside phosphorylase protein (1-289AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Purine nucleoside phosphorylase (PNP) (EC 2.4.2.1) (Inosine phosphorylase) (Inosine-guanosine phosphorylase), PNP, NP
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PNP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.