




TBCEL Antibody

Product Code	CSB-PA722183LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5QJ74
Immunogen	Recombinant Human Tubulin-specific chaperone cofactor E-like protein (1-424AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tubulin-specific chaperone cofactor E-like protein (EL) (Leucine-rich repeat-containing protein 35), TBCEL, LRRC35
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TBCEL

Image



Western blot

All lanes: TBCEL antibody at 12µg/ml

Lane 1: Rat gonadal tissue

Lane 2: PC-3 whole cell lysate

Lane 3: Rat heart tissue

Lane 4: 293T whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 49 kDa

Observed band size: 49 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.