



Serine protease inhibitor Kazal-type 2-like protein Antibody

Product Code	CSB-PA19037A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A0A059T4U4
Immunogen	Recombinant Cyprinus carpio Serine protease inhibitor Kazal-type 2-like protein (18-77AA)
Raised In	Rabbit
Species Reactivity	Cyprinus carpio
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serine protease inhibitor Kazal-type 2-like protein
Immunogen Species	Cyprinus carpio
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.