



# LAG3 Antibody, Biotin conjugated

Product Code	CSB-PA012719LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P18627
Immunogen	Recombinant Human Lymphocyte activation gene 3 protein (29-450AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphocyte activation gene 3 protein (LAG-3) (Protein FDC) (CD antigen CD223), LAG3, FDC
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	LAG3