



Phospholipase D LamSicTox-alphaIC1 Antibody

Product Code	CSB-PA495649LA01LQU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C0JAZ9
Immunogen	Recombinant <i>Loxosceles amazonica</i> Phospholipase D LamSicTox-alphaIC1 protein (1-273AA)
Raised In	Rabbit
Species Reactivity	<i>Loxosceles amazonica</i>
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phospholipase D LamSicTox-alphaIC1 (PLD) (EC 3.1.4.4) (Dermonecrotic toxin) (Sphingomyelin phosphodiesterase D) (SMD) (SMase D) (Sphingomyelinase D) (Fragment)
Immunogen Species	<i>Loxosceles amazonica</i> (Recluse spider)
Research Area	Others
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.