




PTPRZ1 Antibody

Product Code	CSB-PA019068LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P23471
Immunogen	Recombinant Human Receptor-type tyrosine-protein phosphatase zeta protein (36-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor-type tyrosine-protein phosphatase zeta (R-PTP-zeta) (EC 3.1.3.48) (Protein-tyrosine phosphatase receptor type Z polypeptide 1) (Protein-tyrosine phosphatase receptor type Z polypeptide 2) (R-PTP-zeta-2), PTPRZ1, HTPZP2 PTPRZ PTPRZ2 PTPZ
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PTPRZ1

Image



Immunohistochemistry of paraffin-embedded human brain tissue using CSB-PA019068LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.