



PTPRZ1 Antibody, Biotin conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA019068LD01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P23471 |
| Immunogen | Recombinant Human Receptor-type tyrosine-protein phosphatase zeta protein (36-300AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Biotin |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Receptor-type tyrosine-protein phosphatase zeta (R-PTP-zeta) (EC 3.1.3.48) (Protein-tyrosine phosphatase receptor type Z polypeptide 1) (Protein-tyrosine phosphatase receptor type Z polypeptide 2) (R-PTP-zeta-2), PTPRZ1, HTPZP2 PTPRZ PTPRZ2 PTPZ |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Neuroscience |
| Target Names | PTPRZ1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |