



tmk Antibody, HRP conjugated

Product Code	CSB-PA007226HB01BQG
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	C3LJ02
Immunogen	Recombinant Bacillus anthracis Thymidylate kinase protein (1-208AA)
Raised In	Rabbit
Species Reactivity	Bacillus anthracis
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Thymidylate kinase (EC 2.7.4.9) (dTMP kinase), tmk
Immunogen Species	Bacillus anthracis (strain CDC 684 / NRRL 3495)
Research Area	Others
Target Names	tmk
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.